Spectroscopy Spektroskopie

AAS/ICP

## AA-7000F

*Flame system  
Flammensystem*

Three instrument configurations are available, the AA-7000F for analyses in the flame mode, the AA-7000G for the graphite furnace mode and the AA-7000 complete system, featuring both atomization units which are automatically positioned in the optical path when using the AAC-7000 (option).

In combination with the ASC-7000 autosampler, a very high sample throughput can be attained. For the determination of mercury and hydride-forming elements, the HVG-1 is available.

- "high-end" system for flame and graphite furnace modes with automatic atomization unit exchange (requires AA-7000)
- two background compensation methods for solving interference problems in flame- and graphite furnace atomization: Deuterium technique compensates for molecular absorption and particulate caused scattering, high speed self reversal method compensates for spectral interferences by line overlapping
- automatic dilution, matrix modification and addition of buffers

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
*Dual atomizer system / Duales Atomisierungssystem*

## AA-7000

This easy to operate double-beam atomic absorption spectrometer is suitable for fully automatic multi-element analyses and determinations in the ultra-trace range.


The system operates in a wavelength range of 185 - 900 nm and includes a Czerny-Turner monochromator with a holographic grating (1800 lines/mm).

- two background compensation methods for solving interference problems in the entire wavelength range
- multi-element analysis for up to 20 elements in one analysis sequence
- automatic lamp turret with programmable preheating
- software with "Wizard" function for simple operation
- automatic optimization of the heating program
- FDA 21 CFR Part 11 (electronic signature) compliant
- Hardware validation
- Safety functions in flame- and furnace operation

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## Spektroskopie


## ASC-7000

This sample preparation station is suitable for fully automatic operation in the flame- and graphite furnace mode. Eight positions are reserved for standard- and buffer solutions, from which samples can be diluted, mixed and homogenised.

- fully automatic matrix modification in the graphite furnace mode or addition of ionisation buffers in the flame mode (microsampling)
- automatic dilution (ASK-7000)
- 60 samples for 20 elements in the multi-element mode
- no carryover when measuring samples with different concentration

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## HVG-1

The hydride system enables the determination of mercury with the cold-vapour technique as well as the hydride-forming elements (As, Se, Sn, Sb, Te and Bi) with a heated quartz cell.

- sample introduction with 3-channel peristaltic pump
- for continuous measurement
- automation in conjunction with the ASC-7000

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## MVU-1A

The MVU-1A enables determination of mercury with the cold-vapour technique. The sample is reduced in the reaction vessel and the liberated mercury is transferred into the quartz cuvette.

Using this method, a detection limit of 0.01 ppb can be obtained.

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## GFA-7000

The high-sensitivity GFA-7000 graphite furnace enables ultra-trace analysis according to the electrothermal atomization principle. The digital temperature and gas flow control is a new feature which enables optimization of all set parameters for the target element, such as the actual temperature in the drying-, ashing- and atomization stage as well as the heating rate.


- extraordinary graphite tube lifetime for long multi-element sequences
- flexible heating programs with up to 20 temperature steps for complex matrices
- sample enrichment for atomization
- two gas connections for argon or nitrogen and an alternative gas
- furnace camera for optimization of sample injection

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## AA-6200

This cost effective basic instrument for flame emission and atomic absorption is easy to operate. The MS Windows-based user interface offers a vast array of possibilities for system control and data management.

The AA-6200 has double-beam optics for excellent long-term stability, D2-background compensation as well as a 2-lamp turret with lamp preheating.

- compact system with small footprint
- double-beam optics for sensitive measurements
- simple, user-friendly software with "Wizard" function
- electrothermal atomization with GFA-EX7 (optional)
- also optional: ASC-6100 autosampler and GFA-EX7 graphite furnace connection

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## ICPE-9000

ICPE-9000 is a simultaneous ICP-OES spectrometer utilizing a vacuum type Echelle optical system in combination with a newly developed CCD detector. The system is designed for stable performance in the wavelength range from 167 to 800 nm.

In combination with the vertically oriented Mini Torch, argon consumption and therefore the majority of running costs have been significantly reduced.

The ICPE-9000 features easy-to-use, innovative ICPEsolution Software, which offers such features as a Method Development Assistant and a Diagnosis Assistant and easily handles automatic wavelength optimization as well as interference corrections.

- high sample throughput
- low running costs
- stable performance

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**Shimadzu Europa GmbH**