Chromatography Chromatographie

GC/GCMS

GC/GCMS

Shimadzu has been manufacturing gas chromatographs for decades. Customers benefit from this experience by the continuous development of innovative products designed totally for customer needs.

GC-2010 Plus

Separation efficiency, productivity, reduced analysis times – the new GC-2010 Plus offers a new approach in high-end analysis. The system detects and identifies trace-level compounds, supported by state-of-the-art detectors (for instance FID and FPD) featuring world-class sensitivity.

- Trace analysis with high sample throughput enabled by fast GC
- reduced total analysis time due to low inner diameter columns and independent operation of two analytical flow lines

- improved 'double jet' oven cooling system for dynamic heat Abstraction

The new optional AFT (Advanced Flow Technology) functions additionally increase the performance of the GC-2010 Plus. AFT is based on four technologies:

1. MDGC (multidimensional GC) separates and quantifies compounds in difficult matrices
2. a backflush system reduces analysis time by reversing the carrier gas flow and backflushing high volatile compounds

3. a detector splitting system leads the column effluents to multiple detectors, including a mass spectrometer
4. using the detector switching system multiple detectors can be selected easily.

The AFT functions are controlled by the free AFT software.

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Chromatographie

GC-2014

Best results for all budgets and applications – with a range of models for packed column as well as capillary analysis, the GC-2014 is ideally suited for routine analysis.

- electronic control of carrier gas for packed and capillary analysis
- allows Dual flow line systems
- large LCD monitor
- small footprint
- Self-diagnosis functionality and log file; electronic signature with CLASS Agent software (fulfill the requirements of FDA 21 CFR part 11)
- Valves for column switching systems can be mounted

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Shimadzu is also one of the world-wide leading manufacturers of mass spectrometric (MS) equipment and operates a research laboratory for MS headed by Koichi Tanaka who was awarded the Nobel Prize for Chemistry in 2002.

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GCMS-QP2010 Ultra

The newly developed GCMS-QP2010 Ultra surpasses its predecessor model in terms of functionality and user-friendliness. The novel electronics enable up-to-now unmatched data acquisition speed. The system is perfect for fast GC/MS and comprehensive GC/MS (GCxGCxMS). Productivity and sample throughput are increased.

- highest sensitivity measurements made possible by a new ion source and ion optics with highly efficient ion transmission and temperature homogeneity
- increased sensitivity at high-speed data acquisition speeds (up to 20,000 u/s and 100 Hz) via the newly designed and patented ASSP technology (Advanced Scanning Speed Protocol)
- Scan and SIM data within one run enabled by FASST (Fast Automated Scan/SIM Type)
- 'Easy sTop' saves valuable time by enabling injector maintenance without the need to vent the MS
- software-integrated 'Ecology mode' preserves natural resources and decreases costs

- Models featuring EI (electron-impact ionization) as well as models featuring positive and negative chemical ionization (PCI, NCI) are available

Options:

- the Twin Line MS system enables switching between different applications without physically modifying the two installed columns
- Direct-inlet system (DI) for fast information on spectra
- various general MS libraries are available, such as NIST (including AMDIS) and Wiley library as well as specific libraries including retention index information for pesticides (including NCI spectra), fragrances and flavors (FFNSC), FAMEs, drugs and metabolites (MPW and Metabolomics library)
- ChromSquare software for straightforward analysis of comprehensive GCMS measurement data

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GCMS-QP2010 SE


The GCMS-QP2010 SE is a variation of the GCMS-QP2010 Ultra which combines high productivity with exceptional price/performance ratio. It is the ideal tool for routine analysis in laboratories covering different application areas, such as environmental analysis or quality control.

- available as EI model
- Fast GC with new, patented ASSP (Advanced Scanning Speed Protocol) technology with scan speeds up to 10,000 u/s.
- a wide range of accessories supports highly flexible and efficient use. Various general MS libraries are available, such as NIST (including AMDIS) and Wiley library, as well as specific libraries including retention indices for pesticides (including NCI spectra), fragrances and flavors (FFNSC), FAMEs, drugs and metabolites (MPW and Metabolomics library)

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Available Detectors and Injection Systems

FID, (μ -)ECD, FTD (NPD), FPD and (μ -)TCD detectors are available for the GC-2010 Plus and the GC-2014 gas chromatographic systems. All GC-2010 Plus detectors can be installed optionally on the QP2010 series Ultra and SE GCMS systems.

In addition to a split/splitless injector, inlet systems including a combined on-column/PTV (temperature programmable injector), a cooled PTV (Optic 3), as well as a thermodesorption system and a pyrolyzer system are available.

AOC-5000 Plus

The AOC-5000 Plus autosampler for GC and GC/MS provides powerful working capabilities and supports the growth of a laboratory. The system can be upgraded with many modes to meet the challenges of today's requirements.

The AOC-5000 Plus autosampler combines liquid, large volume and headspace injection as well as solid-phase microextraction (SPME) in one single instrument. This unique capability allows quick switching from one application to another.

- superior accuracy and reproducibility combined with optimum injection flexibility
- rapid and efficient sample enrichment of volatile and semivolatile compounds in solid, liquid and gaseous samples
- software control through GCsolution and GCMSsolution using Cycle composer software
- interfaces with the Shimadzu systems GC-2010, GC-2010 Plus, GC-2014, GC-17, GC-14, GCMS-QP2010 series (Plus, Ultra, S, SE) as well as with the GCMS systems QP-5000 and QP-5050
- temperature controlled sample storage from 4 °C up to 40 °C with optional Peltier Stack Cooler
- various options like MHE (Multiple Headspace Extraction), ITEX-2 for efficient extraction of volatile and semi-volatile sample compounds, Dilutor and Flow Cell
- extended version is available

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Chromatography Software

Control of the GC and GCMS systems and data handling are carried out via the easy-to-use GCsolution/GCMSsolution software, part of the LabSolutions series for chromatography. The new ChromSquare software offers straightforward evaluation of comprehensive chromatography data.

AOC-20i/s

High-precision liquid injections – for 6, 12 or 150 samples. The autosampler is self-adjusting and can easily be mounted and is self-adjusting.

- Split/splitless injection, on-column and large-volume injection with multi-injection
- full flexibility for all injection parameters
- all parameters are software controlled (e.g. fill speed, volume, injection time etc.)
- dual autosampler system serves two injectors simultaneously or sequentially

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
TD-20

The new TD-20 thermodesorption system offers optimum performance for the analysis of volatile and semi-volatile compounds. Thermodesorption can be used successfully for many GC and GCMS applications. Typical applications are the determination of organic pollutants in air or the detection of flavor and fragrance compounds in food.

- high recovery rates also for higher boiling point compounds
- minimal carry-over through prevention of cold spots in the system
- excellent reproducibility by electronic flow control (AFC-2010)
- focussing with second trap
- cooling with Peltier element (no coolant required)
- Autosampler with 48 positions
- Control via TD-control software (linked to GC/GCMSsolution software)

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Py-2020iD/iS

Frontier Lab's Py-2020iD/iS pyrolyzer is perfect for the comprehensive analysis of polymer samples, also in routine use with autosampler (up to 48 samples). The double-shot technique enables the analysis of volatile compounds (thermodesorption) and the polymer base (pyrolysis) with the Py-2020iD. Databases for identification are available.

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
Optic 3

The cooled PTV injector (temperature programmable) is the perfect solution for the analysis of volatile compounds, even in trace level concentrations. With DMI (difficult matrix introduction) samples with high matrix content can also be measured. In combination with the AOC-5000 Plus, the LINEX system offers the possibility of an automatic liner exchange on the Optic 3.

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GC/GCMS

To find your local Shimadzu contact
please visit www.shimadzu.eu

Shimadzu Europa GmbH
Albert-Hahn-Str. 6-10 · D-47269 Duisburg

Tel.: +49 - (0)203 - 76 87-0
Fax: +49 - (0)203 - 76 66 25
E-Mail: shimadzu@shimadzu.eu