

P/N 980-04986 Article Structure Description	
Vial type,	1.5ml Screw Neck Vial ND9 Kit includes Vials and MS Caps
Size	32 x 11.6mm
Glass material	Clear glass
Glass quality	1st hydrol. class
Volume	1.5 ml
Cap	Transparent 9mm Short Thread MS Cap with thinned penetration area and diaphragm one component cap without septum, no bleeding, absolutely inert, pierce able like a septa, LC/GC MS certified
Cap center hole size	
Septum type	Septum less seal
Temperature resistant	
Additional remarks	Vial with label and filling lines
Packing	100 vials with caps and septa in PP box

