

P/N 980-08339 Article Structure Description	
Cap	Red PP Short Thread Cap
Cap center hole size	6mm
Septum type	RedRubber/PTFE beige
Septum thickness	1.0 mm
Septum hardness	45° shore A
Temperature resistant	-40 to 110°C
Additional remarks	Cap fitting for Vial 980-08337 RedRubber is actually a silicone rubber. The septum is an alternative to the UltraClean Septum P/N 980-08338 which is recommended for trace, complex, or critical analysis.
Packing	1000 caps and septa in PP box

